



Skipping Devices for Arm Amputees

Whether for fitness or fun, skipping offers a great way to expend energy, and with the use of various attachments allows partial hand and arm amputees to get involved. The following devices make it possible for amputees to enjoy this activity.

Strap device

This skipping rope device made for a partial hand amputee fastens to the wrist and is designed with a custom rope handle.



Spatula

Alternatively, a simple spatula enables partial hand amputees to grab hold of the rope between his/her hand and the pad on the device.

Prosthesis

This custom skipping rope prosthesis is designed for an arm amputee. A device like the one featured here is often designed when an amputee visits a prosthetist and expresses a desire to participate in a certain activity.



TRS Rocky Mountain High Skip Ropes

The **TRS Rocky Mountain High Skip Rope** is intended for use with a prosthesis. The specially designed skip rope adapter is angled to optimize skipping. The adapter screws into a standard prosthetic wrist or disconnect style wrist adapter.

